




EP Consumables Bin and EP Med-waste Bin

Consumables Bin



The EP Consumables Bin places a large volume of materials within easy reach at chairside.

With a variety of adjustable slot-partitions, the Consumables Bin can store up to 10 days of materials required to cover your most common procedures. The unit is easily transported to central Resupply for weekly restocking.

The Bin is available by itself, or paired with our Med-waste bin. Wall and counter brackets are available for either configuration to allow use in any operator.

Standard EP Consumables Bin - Item M1A

- a. Consumables Bin (1)
- b. Aluminum Lid (1)
- c. Flat dividers (6)
- d. Aluminum Wall/Shelf bracket (1)
- e. T-divider (1)
- f. Multi-divider (1)

Med-waste Bin



Med-Waste Bin configuration - Item M1B

- g. Med-waste Bin (1)
- h. Aluminum Lid (1)
- i. Inner removable cover (1)
- j. Sharps Bin (1)

Dimensions	Outside	Inside	Height
Cons. Bin - Top	7.5" x 16.75"	6" x 15"	6.5"
Cons. Bin - Bottom	5.25" x 14.5"		


Med-waste Bin: 10" wide, for an overall Top width of 26.75"



EP Consumables Bin & Med-Waste - Sample Setup



Configured
EP Consumables Bin
& Med-waste with
aluminum wall/shelf
bracket.



- 1 – Cotton Rolls
- 2 – 2x2 Gauze
- 3 – Cotton Tip Applicators
- 4 – Small and Large
Surgical Suction Tips
- 5 – Saliva Ejection,
High-Volume Suction

- 6 – Lower Saliva Ejector
- 7 – Long Disposable Brushes
- 8 – 30-Gauge
Extra-Short Needles
- 9 – 30-Gauge Short Needles
- 10 – Bupivacaine HCl
Anesthetic

- 11 – Articaine HCl Anesthetic
- 12 – 27-Gauge Long Needles
- 13 – Short Disposable Brushes
- 14 – Prilocaine HCl Anesthetic
- 15 – Lidocaine HCl Anesthetic
- 16 – Disposable Mirrors

To see Video of an EP Mega Tub setup, go to <https://youtu.be/iy3Tlv3iBdY>