

## Upgrade Options

*Our delivery platforms offer a wide range of options to suit the way YOU deliver dentistry...*



### Electric Handpiece Integration - Bien Air MX2

We offer full electric integration of the Bien Air MX2 electric handpiece unit. The MX2 brushless micromotor, with adjustable LED light, is simply the best unit available. This revolution in engineering and control technology provides unparalleled balance and comfort, while giving constant power and perfect control of speed and torque.

Control unit and single motor included. Handpieces sold separately.

### Ultrasonic Scaler Integration

Most of our Workstations can incorporate a scaler. We offer a fully integrated Parkell Turbo Sensor, consolidating multiple instrument control in a single foot pedal.



Parkell Turbo Sensor

### Curing Light Options

THE CURE, TC-01, LED light is one of the smallest and most powerful corded lights, providing 1200 mW/cm<sup>2</sup> intensity. The LED has 4 diodes and allows for a 6-second cure.

**The Cure TC-01 by Spring Health** (included as standard on Stratus, as an option on other Workstations)



The CURE, TC-24A, LED light is the latest in high-power LED technology providing 1400 mW/cm<sup>2</sup> intensity. The LED has 4 diodes to provide a 6-second cure. The TC-24A is one of the lightest curing lights available (1.7 oz) and ergonomically designed to provide a comfortable grip from any angle.

**The Cure TC-24A by Spring Health** (available as an upgrade for all Workstations)



### Additional Utility Upgrades Available

- 3rd Air-driven Handpiece, High or Slow-speed (in 4, 5 or 6 hole options)
- 2nd Air/Water Syringe
- Dedicated Air Syringe
- 2nd / 3rd HVE



Not all upgrades are available for all Workstation configurations. See individual item listings for details.